


Zimmer® Contour Zirconia Abutment Restorations



1. During placement, orient the flat of the pre-attached Fixture Mount/Transfer toward the buccal/labial surface. This ensures that the flat of the internal hex is in the correct orientation.




2. Place a healing collar or provisional restoration in a single-stage procedure for soft tissue contouring.



3. Use a light to medium-body material around the transfer and record a full-arch impression with medium-body material. Remove the transfer from the implant after the impression has been taken.




4. Attach an Implant Analog to the transfer with the 1.25mmD (0.050") Hex Tool and re-insert into the impression. Send the impression to the laboratory.



5. Pour the cast in die stone. Use soft tissue material to represent gingival contours. Use the Contour Abutment Try-ins to help select the appropriate angle and cuff height of the final abutment.



6. Mark the Contour Zirconia Abutment for required modifications to achieve appropriate vertical clearance as well as gingival contours.




7. Prepare the zirconia abutment according to the Instructions for Use. Do not reduce the coronal wall of the 3.5mmD zirconia abutment.



8. Modify the abutment with burs specifically designed for use with zirconia. Preparation with water is highly recommended.



9. Prepare the coping using traditional prosthodontic techniques.



10. Complete the final restoration according to traditional laboratory procedures.



11. Place the zirconia abutment and tighten the retaining screw to 30 Ncm with a calibrated prosthetic torque wrench.



12. Block out the screw access channel in the abutment and cement the final prosthesis to complete the restoration.

Note: This guide is for educational use only.

Refer to the Zimmer Dental Instructions for Use for *Zimmer Contour Zirconia Abutment* (part #8905) for complete instructions for use.

