



One Abutment -
One Time™



A new treatment concept –
simple, innovative, easy to implement.

Something so simple that does so much.

What makes *One Abutment - One Time* a favorite among clinicians and patients alike? Simply put, from its ease of implementation to its time-saving convenience, no other abutment solution does so much, for so many, so well.



Benefits

Final abutment is inserted only once.

The friction-fit connection provides a high degree of abutment stability from the onset, virtually eliminating micromovement.¹ The final abutment is inserted immediately, requiring no additional prep work.

Protects soft tissue.

Better for your patients – the immediate restoration protects tissue. Initial gingival attachments to the abutment are no longer destroyed due to repeated placement and removal of abutments and healing abutments, also helping to minimize potential bone loss.


Save time and money.

Placement of conventional healing abutments and subsequent visits for impression-taking are no longer necessary. Prefabricated prosthetic components save time and money in the office and the lab.

Fewer visits for your patients.

Conventional: 4-5 visits

One Abutment - One Time: 2-3 visits
(Visits for placing healing abutments and impression-taking are eliminated.)





Get new referrals.


Simple treatment concept allows the restoring dentist to perform the temporary and final restoration without any additional instruments or investments. Subsequent prosthetic treatment is similar to that of conventional crown-and-bridge techniques.


¹ Binon PP. The evolution and evaluation of two interference-fit implant interfaces. *Postgraduate Dent.* 1996;3:3-13.


One treatment concept, four winners.


- 


1 Initial situation following implant placement and healing phase.
- 

2 Expose and remove the surgical cover screw.
- 

3 Determine exact abutment height using Contour Abutment Try-ins that can be resterilized.
- 

4 Place the *Hex-Lock*® Contour Abutment onto the implant. The short side of the collar should be positioned to the labial or buccal. Use a calibrated torque driver and tighten the abutment screw to 30 Ncm. Take a radiograph to confirm that the abutment is completely seated.
- 

5 Place the Contour Impression Cap over the abutment, making sure the cap is aligned with the contours of the margin. Snap the cap into place.
- 

6 Inject the impression material around the Contour Impression Cap and take a full-arch impression. The cap will be picked up in the impression tray.
- 

7 Prepare the provisional crown by applying acrylic to the Contour Provisional Coping (or use it alone). Block out the screw channel and use temporary cement to cement the crown or Contour Provisional Coping.

One Abutment - One Time.

Dental Lab

1



Align the Contour Abutment Analog with the Contour Impression Cap and place it into the impression, snapping the analog into place.

2



In the lab, pour the impression with die stone. Use a soft tissue material to represent gingival contours.

3



Place the Contour Waxing Coping onto the Contour Abutment Analog in the master cast. The coping does not snap onto the analog.

4



Use traditional prosthetic techniques to wax and cast the coping.

5



Apply ceramic veneering onto the casting to finish the restoration.

6



Fit the final crown onto the master cast.

Patient

1



Remove the provisional crown and any remaining cement. Block out the screw channel and cement the final crown.

2



Final restoration in place.



For more information on our Products and Educational Opportunities, contact us:

1900 Aston Avenue
Carlsbad, CA 92008-7308
USA

In the U.S. 1 (800) 854-7019

To fax an order 1 (888) 225-2483

Outside the U.S. +1 (760) 929-4300

Australia +61 (0)2 9950 5434 or 1 (800) 241 916

Canada + 1 (905) 567-2073 or 1 (800) 265-0968

France +33 (0)1 45 12 35 35

Germany +49 (0)761 1 56 47 0

Israel +972 (0)3 6124242

Italy +39 043 855 5573

Spain +34 93 846 05 43

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsLetterSignUp.aspx



www.zimmerdental.com