




1
Remove healing collars from implants using the 1.25mmD Hex Tool.



2
Screw Tapered Abutments onto implants using the 1.25mmD Hex Tool. Apply 30 Ncm of torque with a calibrated torque wrench. One-piece Tapered Abutments do not engage the internal interface of the implant.



3
Screw Tapered Abutment Indirect Transfers onto abutments using the 1.25mmD Hex Tool. Do not over tighten. An open-tray technique can also be performed using Tapered Abutment Direct Transfers.



4
Use light body impression material around transfers and record a full-arch impression with standard body material.



5
Unscrew transfers from Tapered Abutment and attach the transfers to the abutment replicas. Insert assembly back into impression. If using an open-tray technique, attach replicas to transfers remaining in the impression.



6
Screw healing caps onto abutments using the 1.25mmD Hex Tool.



7
Pour the model in die stone. If desired, use soft tissue material to represent gingival tissues. Remove transfers.



8
Seat Tapered Abutment Gold Copings on to the abutment replicas and secure them with coping screws.



9
Modify sheaths as needed. Wax and cast the framework using traditional prosthodontic techniques. Try the framework on the model.



10
Check framework for accuracy and passive seating in the patient's mouth. If not passive, section and index frame to a passive fit and return to the lab. Try in framework again. Once the passive fit is confirmed, send to laboratory for completion.














11
Remove abutment healing caps. Reapply 30Ncm of torque with a calibrated prosthetic torque wrench. Seat prosthesis and tighten coping screws with 20Ncm torque wrench.



12
Block out screw access holes with material of choice. Final restoration is in place.



Note: This guide is for educational use only. Refer to the Centerpulse Dental Division Prosthetics Instructions for Use (part #4894), and Instructions for Use (part #'s 4869 and 4718) for complete instructions for use.

Tapered Screw-Vent®

Qty.	Tapered Abutments	
	Platform Diameter	Abutment Length
 TAC2	3.5mm	2mm
 TAC3	3.5mm	3mm
 TAC4	3.5mm	4mm
 TAC5	3.5mm	5mm
 TACW2	4.5mm	2mm
 TACW3	4.5mm	3mm
 TACW4	4.5mm	4mm
 TACW5	4.5mm	5mm
 TA5C2	5.7mm	2mm
 TA5C3	5.7mm	3mm
 TA5C4	5.7mm	4mm



AdVent™

Qty.	Tapered Abutments	
	Platform Diameter	Abutment Length
 AVACT	4.5mm	No added height
 AVACT3	4.5mm	3.5mm



Note: Use 5.7mmD components for 5.7mmD platform AdVent.



Tapered Screw-Vent implants feature a patented fixture mount/transfer that can be used for impressions at the time of surgery or prepped as a temporary abutment.

For a complete list of products and tooling, please refer to the current price list.

From the U.S.
 Fax Orders to: (888) 225-2483




From Australia or New Zealand
 Fax Orders to: +61 2 429 74190

From Canada
 Fax Orders to: (416) 674-6122

From Israel
 Fax Orders to: 972 (0)3 6124243











Others
 Fax Orders to: +1(760) 431-7811

Tapered Abutment Transfer Components

Qty.	Description
 ACTIT	Tapered Abutment Indirect Transfer
 ACTDT	Tapered Abutment Direct Transfer
 ACTR	Tapered Abutment Replica





Tapered Abutment Restorative Components™

Qty.	Description
 ACTGC	Tapered Abutment Gold Coping
 ACTP	Plastic Castable Coping
 TGC3	Tapered Abutment Bar Gold Coping - 3mmL
 TGC5	Tapered Abutment Bar Gold Coping - 5mmL
 SCTS	Replacement Coping Screw
 HGB	1.8mm Round Gold Bar with Hader Clips
 DGB	Dolder Resilient Gold Bar
 BS1	Hader Clip Bar System
 SCWS	Tapered Abutment Waxing Screw
 TATHC	Tapered Abutment Healing Cap











SwissPlus™

Qty.	Tapered Abutments	
	Platform Diameter	Abutment Length
 OPACT	4.8mm	1.6mm
 OPACT3	4.8mm	3.0mm






Tools

Qty.	Description
 THX1.25	Tapered Frictional Hex Tool, 1.25mmD
 THXL1.25	Long Tapered Frictional Hex Tool, 1.25mmD
 HX1.25	Hex Tool, 1.25mmD
 HXL1.25	Long Hex Tool, 1.25mmD
 TW30	Prosthetic Torque Wrench, 30Ncm
 TW20	Prosthetic Torque Wrench, 20Ncm
 TW1.25	Torque Wrench Hex Driver, 1.25mmD
 TW1.25L	Long Torque Wrench Hex Driver, 1.25mmD



Other

Qty.	Description
	
	
	

Customer Account _____

Ship To:
 Dr. Name _____

Phone (_____) _____

Address _____

City _____

State _____ Zip _____

Order Date _____

Ordered By _____

Contact Information
 In the U.S. 800 854 7019
 Outside the U.S. + 1 760 929 4300
 France + 33 (0)1 45 12 35 35
 Germany + 49 (0)7 61 45 84 723
 Canada + 1 800 265 0968 or 1 416 674 6116
 Israel + 972 (0)3 6124242
 Australia + 61 2 429 74182 or 1800 241 916

Trademark Information
 Screw-Vent is a registered trademark of Centerpulse Dental Inc.
 Advent and SwissPlus are trademarks of Centerpulse Dental Inc.

Ship Via Federal Express
 2 Day Overnight (5:30 p.m.) P1 (10:30 a.m.)

Must Arrive By _____

Payment Options
 Invoice COD Visa MasterCard AmEx

Credit Card # _____

Exp. Date _____

Signature _____

New Zealand + 61 2 966 20844 or 0800 305 566
 Spain +34 93 846 05 43

www.centerpulse-dental.com

Centerpulse Dental Inc.
 1900 Aston Avenue
 Carlsbad, CA 92008 - 7308
 U.S.A.