



1
In a case with a failing tooth, perform atraumatic extraction, taking care to preserve the alveolar housing around the extraction site.



2
Prepare surgical site according to the drilling sequence to set the angulation of the implant.



3
Deliver the threaded implant to the osteotomy, and seat the implant using a handpiece or manual ratchet.



4
After removal of the mount, select an abutment with the appropriate emergence and angle. Prepare the abutment extraorally and seat the abutment on the implant with a 1.25mmD Hex Tool and a torque wrench calibrated to 30Ncm. The friction-fit abutment can only be attached if the primary stability of the implant is satisfactory. Alternatively prepare the fixture-mount and use as a temporary abutment.



5
Fabricate provisional prosthesis chairside. Check provisional prosthesis to confirm no occlusal contacts are present. Occlude screw access hole and then cement provisional prosthesis in place using temporary cement.



6
After the healing period, remove provisional prosthesis and abutment. Attach a transfer post with a 1.25mmD Hex Tool and record full-arch impression with standard body material. Remove transfer and reseat provisional prosthesis. Forward transfer and impression for fabrication of the final restoration. Optional: Can take a direct impression of abutment.



7
Connect the transfer to a corresponding replica and insert assembly back into the impression. Create working cast using soft tissue material to represent the gingival tissues. Remove transfer from the cast. Attach an abutment to the replica on the cast and tighten the abutment screw to the appropriate torque with a 1.25mmD Hex Tool. Prepare, wax and cast the abutment and fabricate the final prosthesis.



8
Place the prepared final abutment in the patient's mouth. Thread the abutment screw through the access channel within the abutment using the 1.25mmD Hex Tool. Tighten the abutment screw to 30Ncm with a calibrated torque wrench.



9
Occlude screw access hole and then attach the final prosthesis to the prepared abutment with cement of choice.

Note: This guide is for educational use only.
Refer to the Zimmer Dental Division Prosthetics Instructions for Use (part #4894) and Instructions for Use (part #'s 4869 and 4718) for complete instructions for use.

Tapered Screw-Vent® Implants

Qty.	MTX™			
	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	TSVB8	3.7mm	3.5mm	8mm
<input type="radio"/>	TSVB10	3.7mm	3.5mm	10mm
<input type="radio"/>	TSVB13	3.7mm	3.5mm	13mm
<input type="radio"/>	TSVB16	3.7mm	3.5mm	16mm
<input type="radio"/>	TSVWB8	4.7mm	4.5mm	8mm
<input type="radio"/>	TSVWB10	4.7mm	4.5mm	10mm
<input type="radio"/>	TSVWB13	4.7mm	4.5mm	13mm
<input type="radio"/>	TSVWB16	4.7mm	4.5mm	16mm
<input type="radio"/>	TSV6B8	6.0mm	5.7mm	8mm
<input type="radio"/>	TSV6B10	6.0mm	5.7mm	10mm
<input type="radio"/>	TSV6B13	6.0mm	5.7mm	13mm
<input type="radio"/>	TSV6B16	6.0mm	5.7mm	16mm

Qty. MP-1® HA Dual Transition Selective Surface

	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	TSVH8	3.7mm	3.5mm	8mm
<input type="radio"/>	TSVH10	3.7mm	3.5mm	10mm
<input type="radio"/>	TSVH13	3.7mm	3.5mm	13mm
<input type="radio"/>	TSVH16	3.7mm	3.5mm	16mm
<input type="radio"/>	TSVWH8	4.7mm	4.5mm	8mm
<input type="radio"/>	TSVWH10	4.7mm	4.5mm	10mm
<input type="radio"/>	TSVWH13	4.7mm	4.5mm	13mm
<input type="radio"/>	TSVWH16	4.7mm	4.5mm	16mm
<input type="radio"/>	TSV6H8	6.0mm	5.7mm	8mm
<input type="radio"/>	TSV6H10	6.0mm	5.7mm	10mm
<input type="radio"/>	TSV6H13	6.0mm	5.7mm	13mm
<input type="radio"/>	TSV6H16	6.0mm	5.7mm	16mm

Qty. "Cast-To" Gold Abutments, Engaging

<input type="radio"/>	HLA3G	3.5mmD
<input type="radio"/>	HLA4G	4.5mmD
<input type="radio"/>	HLA5G	5.7mmD

Note: 5.7mmD platform prosthetic components are compatible with the *AdVent* 5.7mmD platform.

© 2004 Zimmer Dental Inc. All rights reserved. 4890, Rev. 12/04.

Prosthetics

Qty.	Replicas	
<input type="radio"/>	IA3	3.5mmD
<input type="radio"/>	IA4	4.5mmD
<input type="radio"/>	IA5	5.7mmD

Qty.	Extension Screw for FM/T	
<input type="radio"/>	DHTS	Optional Screw for Open Tray Impression with FMT

Qty.	Hex-Lock™ Abutments		
	Platform Diameter	Profile Diameter	
<input type="radio"/>	HLA3/3	3.5mm	3.5mm
<input type="radio"/>	HLA3/4	3.5mm	4.5mm
<input type="radio"/>	HLA3/5	3.5mm	5.5mm
<input type="radio"/>	HLA4/4	4.5mm	4.5mm
<input type="radio"/>	HLA4/5	4.5mm	5.5mm
<input type="radio"/>	HLA4/6	4.5mm	6.5mm
<input type="radio"/>	HLA5/6	5.7mm	6.5mm

Qty.	20° Angled Abutments for 6 positions		
	Platform Diameter	Profile Diameter	
<input type="radio"/>	AH20/4	3.5mm	4.5mm
<input type="radio"/>	AH20W/5	4.5mm	5.5mm
<input type="radio"/>	A5H20/6	5.7mm	6.5mm

Qty.	AdVent® Prosthetics, 4.5mmD Platform	
<input type="radio"/>	AVIT/4	Indirect Transfer
<input type="radio"/>	AVR	AdVent Replica
<input type="radio"/>	AVHL/4	AdVent Hex-Lock, 4.5mm Flare
<input type="radio"/>	AVHL/6	AdVent Hex-Lock, 6.5mm Flare
<input type="radio"/>	AVH20/4	20° Angled Abutment for 6 positions
<input type="radio"/>	AVGA	"Cast-To" Gold Abutment

AdVent Implants

Qty.	MTX			
	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	AVB8	3.7mm	4.5mm	8mm
<input type="radio"/>	AVB10	3.7mm	4.5mm	10mm
<input type="radio"/>	AVB13	3.7mm	4.5mm	13mm
<input type="radio"/>	AVB16	3.7mm	4.5mm	16mm
<input type="radio"/>	AVWB8	4.7mm	4.5mm	8mm
<input type="radio"/>	AVWB10	4.7mm	4.5mm	10mm
<input type="radio"/>	AVWB13	4.7mm	4.5mm	13mm
<input type="radio"/>	AVWB16	4.7mm	4.5mm	16mm
<input type="radio"/>	AV6B8	4.7mm	5.7mm	8mm
<input type="radio"/>	AV6B10	4.7mm	5.7mm	10mm
<input type="radio"/>	AV6B13	4.7mm	5.7mm	13mm
<input type="radio"/>	AV6B16	4.7mm	5.7mm	16mm

Qty. MP-1 HA Dual Transition Selective Surface

	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	AVH8	3.7mm	4.5mm	8mm
<input type="radio"/>	AVH10	3.7mm	4.5mm	10mm
<input type="radio"/>	AVH13	3.7mm	4.5mm	13mm
<input type="radio"/>	AVH16	3.7mm	4.5mm	16mm
<input type="radio"/>	AVWH8	4.7mm	4.5mm	8mm
<input type="radio"/>	AVWH10	4.7mm	4.5mm	10mm
<input type="radio"/>	AVWH13	4.7mm	4.5mm	13mm
<input type="radio"/>	AVWH16	4.7mm	4.5mm	16mm
<input type="radio"/>	AV6H8	4.7mm	5.7mm	8mm
<input type="radio"/>	AV6H10	4.7mm	5.7mm	10mm
<input type="radio"/>	AV6H13	4.7mm	5.7mm	13mm
<input type="radio"/>	AV6H16	4.7mm	5.7mm	16mm

From the U.S.

Tel: 800 854 7019
Fax Orders to: 888 225 2483

From Australia

Tel: +61 2 4297 4182 or 1 800 241 916
Fax Orders to: +61 2 429 74190

From Canada

Tel: +1 800 265 0968 or 1 905 567 2073
Fax Orders to: 1 905 567 2076

From France

Tel: +33 (0)1 45 12 35 35
Fax Orders to: +33 (0)1 45 60 04 88

From Germany

Tel: +49 (0)761 4584 722/723
Fax Orders to: +49 (0)761 4584 729

From Israel

Tel: +972 (0)3 6124242
Fax Orders to: +972 (0)3 6124243

From Spain

Tel: +34 93 846 05 43
Fax Orders to: +34 93 845 43 25

Others

Tel: +1 760 929 4300
Fax Orders to: +1(760) 431 7811

Customer Account _____

Ship To:

Dr. Name _____

Phone (_____) _____

Address _____

City _____

State _____ Zip _____

Order Date _____

Ordered By _____

Ship Via Federal Express

2 Day Overnight (5:30 p.m.) P1 (10:30 a.m.)

Must Arrive By _____

Payment Options

Invoice COD Visa MasterCard AmEx

Credit Card # _____

Exp. Date _____

Signature _____

