



Zimmer CurV™  
Pre-Shaped  
Collagen  
Membrane

# The latest in RIDGE AUGMENTATION TAKES SHAPE



## 1 Biocompatible

- **Resorbable collagen** does not require removal
- Eliminates invasive, second surgery **compared to non resorbable titanium mesh**

## 2 Focused Bone Augmentation

- Designed to retain graft material during the bone remodeling process, allowing **vertical bone growth** where desired
- Holds particulate to **control and minimize** migration
- **Optimizes ease-of-use and reduces technique-sensitivity**

## 3 Pre-Shaped Configuration

- **Does not require additional structural support**, such as **titanium mesh**
- **Pliable when hydrated**, able to be shaped and trimmed to fit most defect sites
- **Pre-shaped configurations** available for anterior and posterior


# Breaking Free FROM TITANIUM MESH



Zimmer CurV Pre-Shaped Collagen Membrane is resorbable. Designed to retain graft materials such as Puros® Particulate Allograft during ridge augmentation, when fixated with conventional methods (i.e. tacks or screws) - allowing vertical bone growth where desired. Available in both anterior and posterior sizes, it can be trimmed to fit most defect sites.


Because of its inherent resorbability, the Zimmer CurV membrane stands in sharp contrast to titanium mesh. It eliminates the need for an invasive second surgery, and is designed for patient comfort. Zimmer CurV membrane is made from Type 1 collagen derived from bovine Achilles tendon.

## TAKE A CLOSER LOOK (ANTERIOR PLACEMENT)



Clinical photo ©2011 Monish Bhola, DDS, MSD. All rights reserved. Individual results may vary.

## POSTERIOR PLACEMENT




Clinical photos ©2011 Stephen L. Wheeler, DDS. All rights reserved. Individual results may vary.

Membranes are intended to be cut to appropriate size as necessary for the intended surgical application. Measurements are approximate.

## Ordering information

Catalog Number	Description
0174	Posterior
0175	Anterior



To learn more about the Zimmer CurV membrane, please visit us online at [www.zimmerdental.com/nomoremesh](http://www.zimmerdental.com/nomoremesh), or to speak to a sales representative, call 1 (800) 854 7019.

To receive our eNews visit us at [http://www.zimmerdental.com/news\\_eNewsletterSignUp.aspx](http://www.zimmerdental.com/news_eNewsletterSignUp.aspx)

For more information about our Products, Regenerative Materials and Educational Opportunities, contact us:

In the U.S. 1 (800) 854-7019  
To fax an order 1 (888) 225-2483  
Outside the U.S. +1 (760) 929-4300  
Australia +61 (0)2 9950 5434 or 1 (800) 241 916  
Canada + 1 (905) 567-2073 or 1 (800) 265-0968  
China + 86 21 22115147

France +33 (0)1 45 12 35 35  
Germany +49 (0)761 1 56 47 0  
Israel +972 (0)3 6124242  
Italy +39 043 855 5573  
Spain +34 93 846 05 43



Zimmer Dental  
1900 Aston Avenue  
Carlsbad, CA 92008-7308  
USA

[www.zimmerdental.com](http://www.zimmerdental.com)

©2011 Zimmer Dental Inc. All rights reserved. 2571, Rev. 7/11. Zimmer CurV membrane is manufactured by Osseous Technologies of America for Zimmer Dental Inc. Puros Allografts are manufactured by RTI Biologics, Inc. and marketed by Zimmer Dental Inc. Please note that not all products and regenerative materials are registered or available in every country/region, and part numbers for Puros Allografts vary. Please check with a Zimmer Dental representative for availability and additional information.