

Dependable precision. Enduring performance.



**Zimmer Dental**  
1900 Aston Avenue  
Carlsbad, CA 92008-7308  
USA

[www.zimmerdental.com](http://www.zimmerdental.com)

To receive our eNews visit us at [http://www.zimmerdental.com/news\\_eNewsLetterSignUp.aspx](http://www.zimmerdental.com/news_eNewsLetterSignUp.aspx)

For more information about our Products, Professional Programs  
and Continuing Education, contact us:

In the U.S. 800 854 7019  
To fax an order 888 225 2483  
Outside the U.S. +1 760 929 4300  
Australia +61 2 4297 4182 or 1800 241 916  
Canada +1 800 265 0968 or +1 905 567 2073

France +33 (0)1 45 12 35 35  
Germany +49 (0)761 4584 722/723  
Israel +972 (0)3 6124242  
Spain +34 93 846 05 43

## Starts sharp. Stays sharp.

Introducing new *Dríva* reusable drills from Zimmer Dental, a world leader in dental implant systems. Extensive research and testing have produced a superior drill offering greatly improved sharpness and longevity. Using an advanced grade of surgical stainless steel, *Dríva* drills offer unrivaled levels of strength and cutting efficiency for a more precise and accurate procedure for the life of the drill. The progressive design features a patent-pending cutting geometry for enhanced precision and cutting ability, along with a state of the art coating that reduces glare and improves visualization. The result is simply revolutionary.



## Innovation. Strength. Longevity.

High-performance *Dríva* drills offer clear-cut advantages.

### Advanced cutting efficiency


A sharper, stronger drill. Advanced grade of stainless steel and progressive design with improved cutting geometry makes osteotomy formation faster and more precise.

### Improved longevity

Extensive research and testing for improved endurance and performance. Revolutionary corrosion resistant coating protects drill sharpness and increases efficacy over time. A longer lasting drill saves money.

### Enhanced visibility

Laser-etched depth bands facilitate visual measurement along with a matte finish that reduces glare for an optimal surgical experience.



- Insertion force required with custom stainless steel and advanced coated drill
- Insertion force required with traditional stainless steel and non-coated drill

During in-vitro simulation over twenty-five surgical sites per drill, *Dríva* drills using advanced surgical materials required an average of 41% less insertion force to create an osteotomy when measured against comparable drills manufactured from traditional surgical materials.

Proudly offering one of the most comprehensive dental implant product lines available, Zimmer Dental is a market leader in the development of world-class implantology products, practice-building strategies and educational programs focused on empowering clinicians and improving patients' lives.

*Dríva* drills are available for *Tapered Screw-Vent*®, *Screw-Vent*®, *AdVent*®, *Spline*®, and *SwissPlus*® implants. Contact your Zimmer Dental Sales Consultant or Customer Service at 800 854 7019 for more information, or to place your order.