




## Simple and esthetic provisional restorations



Implanting trust.

# Hex-Lock Plastic Temporary Abutment

Catalog Numbers and Pricing						
		Description	Emergence Profile	Catalog Number	Price	
			3.5mm Hex-Lock Plastic Temporary Abutment	4.5mm	HLPT3	\$40
			4.5mm Hex-Lock Plastic Temporary Abutment	5.5mm	HLPT4	\$40
			5.7mm Hex-Lock Plastic Temporary Abutment	6.5mm	HLPT5	\$40



Hex-Lock Plastic Temporary Abutment includes a Retaining Screw.

- Create a natural emergence profile during tissue healing with simple and esthetic provisional restorations.
- Cement a restoration onto the abutment with temporary cement, or bond autopolymerizing acrylic directly to the abutment.
- A tooth-colored shade enhances esthetics, ideal for anterior restorations.
- Packaged with an easily prepped Retaining Screw that is useful to maintain the screw access channel.
- Affordable for routine use.

For more information about our Products, Professional Programs and Continuing Education, contact us:

In the U.S. 800 854 7019  
 To fax an order 888 225 2483  
 Outside the U.S. +1 760 929 4300  
 Germany +49 (0)7 61 45 84 723  
 Australia +61 2 966 20844 or  
 1800 241 916  
 France +33 (0)1 45 12 35 35  
 Canada +1 800 265 0968 or  
 1 416 674 6116  
 Israel +972 (0)3 6124242  
 New Zealand +61 2 966 20844 or  
 0800 305 566

[www.centerpulse-dental.com](http://www.centerpulse-dental.com)

Centerpulse Dental Inc.  
 1900 Aston Avenue  
 Carlsbad, California 92008-7308  
 U.S.A.